


RF+EMS EYE MASSAGER



Charging Voltage > 5V
Rf frequency > 150kHz
Charging time > 2H
Rated Current > 500mA

✉ hello@bolva-co.com www.bolva-co.com

● Platinum color

● gradient powder



RF+EMS EYE MASSAGER

RF

Our innovative Radio Frequency (RF) Skin Rejuvenation technology delivers 6X more skin tightening with deep energy wave penetration, leaving skin looking younger and firmer than ever before.

EMS

EMS mode harnesses the power of microcurrent technology to lift, tighten and tone skin in minutes. This innovative facial toning system stimulates muscles and tissue with gentle electronic waves to beautify skin from the inside out.



THREE LEVELS



Tightening

massaging for eye muscles to relieve eye fatigue and reduce puffiness. Significantly improve absorption.



Rejuvenating

Stimulating collagen reproduction to erase fine lines and reduce dark circles.



Brightening

Deep tissue lifting to tighten and tone the appearance of skin and remove eyebags.



High Frequency
150kHz



24K Gold
Massage Probes



3 Seconds Warming
and 42°C Temperature Control



Built-in LED
Light Therapy